

Original strap winders since 2014

[dynello.com](http://dynello.com)



**OTTENYD**<sup>®</sup>


**DYNELLO**<sup>®</sup>

**DYNELLO<sup>®</sup> Bandage Roller Drill Bit**

DYNELLO<sup>®</sup> Bandagenwickler Bit

Embout/Bit Enrouleur de Bandage DYNELLO<sup>®</sup>

**RS0305**



# INSTRUCTION MANUAL

BEDIENUNGSANLEITUNG, MANUEL D'INSTRUCTIONS



## Intended use

This **attachment** has been designed for winding applications of lightweight bandages.

## Applicable power units

Only use the attachment along with a power unit with a maximum rotational speed setting of 250 RPM.

## Recommended applicable power unit

It is recommended to use the attachment with a **cordless drill type 3.6V** with maximum rotational speed of 250 RPM.

## Operating tips

Use your free hand to steer the bandage as it is being rolled up, to ensure that it wraps tightly. **Caution!** Do not hold or pull the bandage while it is being rolled up, as this could damage the attachment.

